

Fibrinogen Standard

REF 354 001 1 X 0.5 ml

SYMBOLS IN PRODUCT LABELLING			
	Authorised Representative		Temperature Limitation
	For in-vitro diagnostic use		Use by/Expiration Date
	Batch Code/Lot number		CAUTION. Consult instructions for use
	Catalogue Number		Manufactured by
	Consult instructions for use		

Intended Use

Preparation of reference curves for quantitative immunochemical determination of Fibrinogen in citrated plasma by turbidimetry and nephelometry.

composition

A lyophilized preparation of stabilized, buffered human plasma containing less than 5% bovine plasma.

Preparation

1. Reconstitute the Fibrinogen standard with 0.5 ml of distilled or deionised water.
2. Restopper vial and allow to stand until fully hydrated (a minimum of 5 minutes).
3. Gently mix by rotation or inversion until dissolution is complete.
DO NOT SHAKE.
4. Allow to equilibrate for a total of 30 minutes before use.

Storage and Stability

Stored at +2 to +8 C the Fibrinogen standard remains stable up to labelled expiry date. After reconstitution, the Fibrinogen standard is stable for 48 hours at +2 to +8 C when stored in the original vial in a sealed polypropylene container.

Precautions and Warnings

1. Fibrinogen standard is for "In Vitro Diagnostic Use"
2. Fibrinogen standard is a potentially biohazardous material.
Source materials from which this product was derived were found negative for HBsAg and for antibodies against HCV, HIV-1 and HIV-2 by approved test methods. However, since no test method can offer complete assurance that infectious agents are absent, this product should be handled observing the same safety precautions employed when handling any potentially infectious material.
3. Fibrinogen Standard is Harmful.
Harmful in contact with skin. Irritating to eye, respiratory system and skin. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

Assigned Value

The value is lot dependent. Actual value is printed on the label.



Egyptian Company for Biotechnology (S.A.E)
 Obour city industrial area. block 20008 piece 19 A. Cairo. Egypt.
 Tel: +202 4665 1848 - Fax: +202 4665 1847
 www.spectrum-diagnostics.com
 E-mail: info@spectrum-diagnostics.com



MDSS GmbH
 Schiffgraben 41
 30175 Hannover, Germany



IFUFTI34

Rev.(2), 15/1/2015

