

RHEUMATOID FACTOR (RF) Turbi Latex

REF: 561 001 100 test
R1 Buffer 2 x 20 ml
R2 Latex 1 x 10 ml

Intended Use

In vitro diagnostic reagents for the quantitative determination of Rheumatoid Factor (RF) in human serum by means of particle-enhanced turbidimetric immunoassay.

Background

The most consistent serological feature of rheumatoid arthritis is the increased concentration of autoantibodies directed against antigenic sites in the Fc region of human and animal IgG, namely rheumatoid factors (RFs) in the blood and joint fluid. The potential role of these factors in the pathogenesis of this disease has been studied extensively, with the finding that both environmental and genetic factors affect production of RF. RF determinations are clinically important for the diagnosis, prognosis, and assessment of therapeutic efficacy of rheumatoid arthritis. Although RFs may be found in all immunoglobulin classes, the RF most frequently detected in the laboratory is IgM type, present in about 75 - 80 % of adult patients with rheumatoid arthritis but in about 10 % of children with juvenile rheumatoid arthritis.

Test Principle

This RF test is based upon the reactions between IgM-anti-IgG (RF) in patient's sample and latexcovalently bound human IgG. RF values are determined photometrically.

Reagents

Buffer

Phosphate buffer (0,05 M) pH: 8,0 containing NaCl (0,15M), detergent and polyethyleneglycol.

Preservative : sodium azide < 1g/L

Latex reagent

A suspension of latex microparticules covalently bound human IgG in a glycine buffer (0,1 M, pH: 8,2), containing NaCl (0,15 M) and bovine serum albumin (0,5%).

Preservative: Sodium azide 0,075%

Calibrator

Human - based reference fluid. Preservative: sodium azide, 0.075 %.

All raw materials of human origin used in the manufacture of this product showed no reactivity when tested for HBsAg, anti-HIV-1/2 and HCV with commercially available test methods. However, this product should be handled as though capable of transmitting infectious diseases

preparation




RF Calibrator : Reconstitute with 2 ml of distilled water : Mix gently and bring to room temperature for about 10 minutes before use

Calibration Curve :

prepare the following RF Calibrator dilutions in NaCl 9 g/L. Multiply the concentration of the RF Calibrator by the corresponding factor stated in table below to obtain the RF concentration of each dilution.

Calibrator dilution	1	2	3	4	5	6
Calibrator RF	---	25	50	100	200	400
Na Cl 9 g/L	400	375	350	300	200	---
Factor	0	0.0625	0.125	0.25	0.5	1
Concentration(IU/ml) (for example: the undiluted C = 200 IU/ml)	0	12.5	25	50	100	200

SYMBOLS IN PRODUCT LABELLING

	Authorised Representative		Temperature Limitation
	For in-vitro diagnostic use		Use by/Expiration Date
	Batch Code/Lot number		CAUTION. Consult instructions for use
	Catalogue Number		Manufactured by
	Consult instructions for use		

Storage and Stability

All the components of the kit are stable until the expiration date on the label when stored tightly closed at (2 - 8 °C) and contaminations prevented during their use .Do not use reagents over the expiration date .

Reagent deterioration : presence of particles and turbidity
Reconstituted calibrator :stable for 1 month at 2 - 8 °C or 3 months at 20 °C

Do not freeze: frozen Latex or Diluent could change the functionality of the test.

Samples

Fresh serum or plasma . stable 7 days at2 - 8 °C or 3 months at 20 °C the samples with presence of fibrin should be centrifuged before testing Do not use highly hemolized or lipemic samples .

Procedure

1 - Bring the reagents and the photometer to 37°C

2 - Assay conditions:

Wavelength 650 nm (600 -650 nm)

Temperature 37°C

Cuvette 1cm light path

3 - Adjust the instrument to zero with distilled water .

4 - Pipette into a cuvette :

	blank
Diluent (ml)	0.4
Latex (ml)	0.1

5 - Mix and read the absorbance (blank reagent)

6 - add the sample / Calibrator

	Blank	Sample / Calibrator
NaCl 9 g/L	4 µl	--
Sample / Calibrator	--	4 µl

7 - Mix and read the absorbance after 2 minutes (A2) of the sample addition.

Calculation

Calculate the absorbance difference

(A2-A blank) of each point of the calibrator dilution and plot the values obtained against the RF concentration of each calibrator dilution . Rheumatoid factor concentration in the sample is calculated by interpolation of its (A2-Abblank) in the calibration curve.

Quality Control

Control sera are recommended to monitor the performance of manual and automated assay procedures .

Each laboratory should establish its own Quality Control Scheme and corrective actions if controls do not meet the acceptable tolerances

Expected Values

Up to 20 IU/mL.

Each laboratory should establish its own reference range.

Sensitivity

6 IU /mL.

Linearity

160 IU /mL.
specimens showing higher concentration should be diluted 1+4 using physiological saline and repeat the assay (result×5).

interferences

Hemoglobin (10 g/L) , bilirubin (20 mg/dL) and lipemia (10 g/L) , do not interfere. Other substances may interfere.

References

Arnet FC, Edworthy SM, Bloch DA, McShane DJ, et al. The American Rheumatism Association 1987. Revised criteria for the classification of rheumatoid arthritis. Arthritis Rheum 1988; 31:315-24.

Bartfield H. Distribution of rheumatoid factor activity in non rheumatoid states. Ann NY Acad Sci 1969; 168:30-40. Singer JM, Plotz CM. The latex fixation test. Am J Med 1956; 21:888-92.

Moore TL, Dorner RN. Rheumatoid factors. Clin Biochem 1993; 26:75-84.

Young DS. Effects of Drugs on Clinical Laboratory Test. 5th Edition, AACC Press, 2000.

Passing H, Bablok W. A new biometrical procedure for testing the equality of measurements from two analytical methods.

Application of linear regression procedures for method comparison studies. Part I. J Clin Chem Clin Biochem 1983; 21:709-20.

Sonderdruck aus DG Klinische Chemie Mitteilungen 1995; 26: 207 – 224

ORDERING INFORMATION	
CATALOG NO.	QUANTITY
561 001	100 test



Egyptian Company for Biotechnology (S.A.E)

Obour city industrial area. block 20008 piece 19 A. Cairo. Egypt.
Tel: +202 4665 1848 - Fax: +202 4665 1847
www.spectrum-diagnostics.com
E-mail: info@spectrum-diagnostics.com



MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany



IFUFT124

Rev.(4), 28/4/2014