

## Hemoglobin A<sub>1c</sub> (HbA<sub>1c</sub>) Turbidimetric Immunoassay

REF: 602 001     50 test  
Reagent1     1 x 20 ml  
Reagent2     1 x 4 ml

### Intended Use

Spectrum Diagnostics Hemoglobin A<sub>1c</sub> reagent is intended for Quantitative turbidimetric determination of HbA<sub>1c</sub> in human blood.

### Background

The glycemic control in diabetes mellitus is mainly by the determination of glucose, but also through quantitative determination of hemoglobin A<sub>1c</sub> in human blood. HbA<sub>1c</sub> is an indication for the actual glucose levels over the preceding 3 months. It was shown that HbA<sub>1c</sub> in diabetic subjects can be elevated 2-3 fold over normal and on other hand approaches normal values when they are under metabolic control.

### Assay Principle

This method utilizes the interaction of antigen and antibody to determine the HbA<sub>1c</sub> in whole EDTA blood. HbA<sub>1c</sub> in test samples is adsorbed onto the surface of latex particles, which react with Anti-HbA<sub>1c</sub> (antigen-antibody reaction) and gives agglutination. The amount of agglutination is measured as absorbance. The HbA<sub>1c</sub> value is obtained from a calibration curve.

### Reagent

#### Reagent1 (R1)

Latex.  
Sodium azide (0.95 g/L).

#### Reagent2 (R2)

Anti-human hemoglobin A<sub>1c</sub> mouse monoclonal antibody.  
Stabilizers.

### Materials required but not provided with the kit

#### 1- Standard set

HbA<sub>1c</sub> concentration is stated on the vials labels.

#### 2-Controls

### Reagent Preparation, Storage and Stability

Spectrum HbA<sub>1c</sub> reagents are stable up to the expiry date labeled on the bottles when stored at 2 - 8°C and contaminations are prevented during their use. Once opened the reagents are stable for 1 month if stored tightly closed at 2 - 8 °C after use.

### Specimen Collection and Preparation

Fresh EDTA blood.

### Hemolysate procedure

To determine HbA<sub>1c</sub>, a hemolysate must be prepared for each sample as follow:










1. Dispense 2 ml hemolysis reagent into a test tube.
2. Place 20 µl of well mixed whole EDTA blood (Samples, Standards and Controls) into the test tube and mix.
3. Allow to rest 5 minutes or until complete lysis is evident. Stability of the hemolysate: 72 hours at 2 - 8°C.

### Procedure

Wavelength     650 nm  
Temperature     37 °C  
Cuvette     1cm light path  
Zero adjustment     distilled water

### Solve and lyse standard/control

### SYMBOLS IN PRODUCT LABELLING

	Authorised Representative		Temperature Limitation
	For in-vitro diagnostic use		Use by/Expiration Date
	Batch Code/Lot number		CAUTION. Consult instructions for use
	Catalogue Number		Manufactured by
	Consult instructions for use		

	Standard	Sample
Reagent (R1)	375 µl	375 µl
Standard	5 µl	-----
Sample	-----	5 µl

Mix, and incubate for 2 minutes, then add

Reagent (R2)	75 µl	75 µl
--------------	-------	-------

Mix and read absorbance (A1), incubate for 5 minutes and read absorbance (A2)

Adaptation sheets for several automatic analyzers are available upon request.

### Calculation

Generate a reference curve using HbA<sub>1c</sub> standard set. Determine Δ absorbance of the sample and each standard as following:  
Δ absorbance of sample = (A2 - A1) sample  
Δ absorbance of each standard = (A2 - A1) for each Standard  
Plot the calibration curve and obtain the result.

### Expected Values

Non-diabetics     < 6 %  
Therapeutic diabetics     < 7 %  
Each laboratory should establish its own reference range.

### Linearity

Up to 15 %.  
specimens showing higher concentration should be diluted 1/5 using physiological saline and repeat the assay.

### Dynamic Range

0 - 15 %.

### Waste Disposal

This product is made to be used in professional laboratories. Please consult local regulations for a correct waste disposal.

**S56:** dispose of this material and its container at hazardous or special waste collection point.

**S57:** use appropriate container to avoid environmental contamination.


**S61:** avoid release in environment. refer to special instructions/safety data sheets.

### References

1. Bates, H.M., Lab. Mang., Vol 16 (Jan. 1978)
2. Gonen, B., and Rubenstein, A.H., Diabetologia 15, 1 (1978).
3. Trivelli, L.A., Ranney, H.M., and Lai, H.T., New eng. J. Med. 284, 353 (1971).

### ORDERING INFORMATION

CATALOG NO.	QUANTITY
602 001	50 test

 **Egyptian Company for Biotechnology (S.A.E)**  
Obour city industrial area. block 20008 piece 19 A. Cairo. Egypt.  
Tel: +202 4665 1848 - Fax: +202 4665 1847  
www.spectrum-diagnostics.com  
E-mail: info@spectrum-diagnostics.com

 **MDSS GmbH**  
Schiffgraben 41  
30175 Hannover, Germany



IFUFTI37

Rev.(2), 30/11/2014