

# TWIN A1C

**Ultimate Formula**  
to Monitor Diabetes



**HbA1c and eAG** in one sample

**Intelligent calibration using radio frequency technology**

- No interference from Hb Variants
- Cost effectiveness per test
- Easy access to PC
- Patient archiving software
- Accurate result in 7 minutes
- Portable with built-in battery




	HbA1c Test Score	Estimated Average Glucose	
		mg/dl	mmol/L
Action Suggested	14.0	380	21.1
	13.0	350	19.3
	12.0	315	17.4
	11.0	280	15.6
	10.0	250	13.7
Good	8.0	180	10.0
	7.0	150	8.2
Excellent	6.0	115	6.3
	5.0	80	4.7
	4.0	50	2.6

Diabetes Control Charts


## TWIN A1c Test



Blood Sample Preparation



Hemolysate Preparation



Mix Carefully



Locate in reading location & press enter




After 2 min add 75 µL of R2



Record the result

Easy  
access  
to PC



**Egyptian Co. for Biotechnology (Spectrum Diagnostics)**

Obour city industrial area. Block 20008, piece 19 A. PO Box 30 Obour city

Tel: +202 44892248 / Fax: +202 44892247

Email: info@spectrum-Diagnostics.com - Web: www.spectrum-Diagnostics.com

